


Phospho-c-Kit Antibody (pY362) (F48716)

Catalog No.	Formulation	Size
F48716-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F48716-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P10721
Applications	Dot Blot : 1:500
Limitations	This phospho-c-Kit antibody is available for research use only.

Dot blot analysis of phospho-c-Kit antibody. 50ng of phos-peptide or nonphos-peptide per dot were spotted.



Description

This gene encodes the human homolog of the proto-oncogene c-kit. C-kit was first identified as the cellular homolog of the feline sarcoma viral oncogene v-kit. This protein is a type 3 transmembrane receptor for MGF (mast cell growth factor, also known as stem cell factor). Mutations in this gene are associated with gastrointestinal stromal tumors, mast cell disease, acute myelogenous leukemia, and piebaldism. Multiple transcript variants encoding different isoforms have been found for this gene.

Application Notes

Titration of the phospho-c-Kit antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This phospho-c-Kit antibody was produced from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding pY362 of human c-Kit.

Storage

Aliquot the phospho-c-Kit antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.