


# A4GNT Antibody (R35144)

Catalog No.	Formulation	Size
R35144-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51146
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This A4GNT antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Alpha-1,4-N-acetylglucosaminyltransferase

### **Application Notes**

Optimal dilution of the A4GNT antibody should be determined by the researcher.

#### **Immunogen**

Amino acids PRTYRDLIKGPEGS were used as the immunogen for this A4GNT antibody. **Storage** Aliquot and store the A4GNT antibody at -20oC.