


AAAS Antibody / Adracalin (F53900)

Catalog No.	Formulation	Size
F53900-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53900-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9NRG9
Applications	Western Blot : 1:1000-2000
Limitations	This AAAS antibody is available for research use only.



Western blot testing of human 1) 293T/17 and 2) HeLa cell lysate with AAAS antibody at 1:2000. Predicted molecular weight: 60 kDa.

Description

Achalasia adrenocortical insufficiency alacrimia plays a role in the normal development of the peripheral and central nervous system.

Application Notes

Titration of the AAAS antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 395-425 from the human protein was used as the immunogen for the AAAS antibody.

Storage

Aliquot the AAAS antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.