


# AARS2 Antibody (R35297)

Catalog No.	Formulation	Size
R35297-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57505
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AARS2 antibody is available for research use only.



Western blot of 143B cells overexpressing human AARS2 tested with AARS2 antibody (mock transfection in lane 2). Predicted molecular weight: ~107kDa.

# **Description**

Additional name(s) for this target protein: Alanyl-tRNA synthetase 2

## **Application Notes**

Optimal dilution of the AARS2 antibody should be determined by the researcher.

### **Immunogen**

Amino acids QHRARQAEPVQK were used as the immunogen for this AARS2 antibody.

Storage Aliquot and store the AARS2 antibody at -20oC.						