


# ABHD3 Antibody (F43291)

Catalog No.	Formulation	Size
F43291-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43291-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8WU67
Applications	Western Blot : 1:1000
Limitations	This ABHD3 antibody is available for research use only.



#### **Description**

Abhydrolase domain-containing protein 3 contains an alpha/beta hydrolase fold, which is a catalytic domain found in a very wide range of enzymes. The function of this protein has not been determined.

#### **Application Notes**

Titration of the ABHD3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### **Immunogen**

A portion of amino acids 282-310 from the human protein was used as the immunogen for this ABHD3 antibody.

## **Storage**

Aliquot the ABHD3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.