

## ABTB2 Antibody (F42258)

Catalog No.	Formulation	Size
F42258-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F42258-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q8N961
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This ABTB2 antibody is available for research use only.



ABTB2 antibody western blot analysis in human K562 lysate. Predicted molecular weight ~114/93 kDa (isoforms 1/2).

## Description

Ankyrin repeat and BTB/POZ domain-containing protein 2 may be involved in the initiation of hepatocyte growth (By similarity).

## Application Notes

Titration of the ABTB2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 681-710 from the human protein was used as the immunogen for this ABTB2 antibody.

## Storage

Aliquot the ABTB2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.