


ADAMTS17 Antibody (F54401)


Catalog No.	Formulation	Size
F54401-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54401-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity purified
UniProt	Q8TE56
Localization	Nuclear
Applications	Immunohistochemistry (FFPE) : 1:25 Western Blot : 1:500-1:2000
Limitations	This ADAMTS17 antibody is available for research use only.



IHC testing of FFPE human lung tissue with ADAMTS17 antibody. HIER: steam section in pH9 EDTA for 20 min and allow to cool prior to staining.



Western blot testing of human 1) A431, 2) HEK293, 3) SK-BR-3, 4) RPMI-8226, 5) MOLT4 and 6) PC-3 cell lysate with ADAMTS17 antibody. Predicted molecular weight ~121 kDa.

Description

Alternate name(s): A disintegrin and metalloproteinase with thrombospondin motifs 17

Application Notes

The stated application concentrations are suggested starting points. Titration of the ADAMTS17 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 508-540 from the human protein was used as the immunogen for the ADAMTS17 antibody.

Storage

Aliquot the ADAMTS17 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.