


AHSG Antibody (R35387)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35387-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|--------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 197 |
| Applications | Western Blot : 0.2-0.5ug/ml ELISA (peptide) LOD : 1:128000 |
| Limitations | This AHSG antibody is available for research use only. |



Western blot testing of human albumin-depleted serum with AHSG antibody at 0.5ug/ml.
Predicted molecular weight: ~39/45-59kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Alpha-2-HS-glycoprotein, fetuin-A

Application Notes

Optimal dilution of the AHSG antibody should be determined by the researcher.

Immunogen

Amino acids QPEGANEAVPTP were used as the immunogen for this AHSG antibody.

Storage

Aliquot and store the AHSG antibody at -20oC.