


ALDOA Antibody (R35313)

Catalog No.	Formulation	Size
R35313-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	226
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ALDOA antibody is available for research use only.



Western blot testing of A) human, B) mouse and C) rat skeletal muscle lysates with Aldolase A antibody at 0.03ug/ml. Predicted molecular weight: ~40kDa.

Description

Additional name(s) for this target protein: Aldolase A

Application Notes

Optimal dilution of the ALDOA antibody should be determined by the researcher.

Immunogen

Amino acids QKADDGRPFPQ were used as the immunogen for this ALDOA antibody. **Storage** Aliquot and store the ALDOA antibody at -20oC.