


ANGPTL6 Antibody (F50753)

| Catalog No. | Formulation | Size |
|---------------|--|---------|
| F50753-0.4ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.4 ml |
| F50753-0.08ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.08 ml |

Bulk quote request

| Availability | 1-3 business days |
|--------------------|---|
| Species Reactivity | Human |
| Format | Purified |
| Clonality | Polyclonal (rabbit origin) |
| Isotype | Rabbit Ig |
| Purity | Purified |
| UniProt | Q8NI99 |
| Applications | Western Blot: 1:1000 |
| Limitations | This ANGPTL6 antibody is available for research use only. |



Description

ANGPTL6 may play a role in the wound healing process. It may promote epidermal proliferation, remodeling and regeneration and may promote the chemotactic activity of endothelial cells and induce neovascularization. ANGPTL6 may counteract high-fat diet-induced obesity and related insulin resistance through increased energy expenditure.

Application Notes

Titration of the ANGPTL6 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 206-236 from the human protein was used as the immunogen for this ANGPTL6 antibody. **Storage** Aliquot the ANGPTL6 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com Copyright $\ensuremath{\texttt{@}}$ NSJ Bioreagents. All rights reserved.