

APOA2 Antibody (F52381)

Catalog No.	Formulation	Size
F52381-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52381-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Primate
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P02652
Applications	Western Blot : 1:1000
Limitations	This APOA2 antibody is available for research use only.



APOA2 antibody western blot analysis in human blood plasma lysate. Predicted molecular weight ~11 kDa.

Description

May stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism.

Application Notes

Titration of the APOA2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 23-56 from the human protein was used as the immunogen for this APOA2 antibody.

Storage

Aliquot the APOA2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.