


APOBEC2 Antibody (R33695)

Catalog No.	Formulation	Size
R33695-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10930
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This APOBEC2 antibody is available for research use only.



Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 2

Application Notes

Optimal dilution of the APOBEC2 antibody should be determined by the researcher.

Immunogen

Amino acids EAASQNGEDLENLDD were used as the immunogen for this APOBEC2 antibody. **Storage** Aliquot and store the APOBEC2 antibody at -20oC.