


# **APOBEC3D Antibody (R35050)**

Catalog No.	Formulation	Size
R35050-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	140564
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This APOBEC3D antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme catalytic polypeptide-like 3D

### **Application Notes**

Optimal dilution of the APOBEC3D antibody should be determined by the researcher.

### **Immunogen**

Amino acids PKRQSNHRQEVYFR were used as the immunogen for this APOBEC3D antibody.

# **Storage** Aliquot and store the APOBEC3D antibody at -20oC.