


# APOL2 Antibody (R36192)

Catalog No.	Formulation	Size
R36192-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23780
Applications	Western Blot: 0.03-0.1ug/ml ELISA (peptide) LOD: 1:64000
Limitations	This APOL2 antibody is available for research use only.



# **Description**

Additional name(s) for this target protein: Apolipoprotein L 2; APOL-II

### **Application Notes**

Optimal dilution of the APOL2 antibody should be determined by the researcher.

## **Immunogen**

Amino acids TKIHEMLQPGQD were used as the immunogen for this APOL2 antibody.

Storage Aliquot and store the APOL2 antibody at -20oC.				