


## APP Antibody (R33373)

Catalog No.	Formulation	Size
R33373-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	351
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This APP antibody is available for research use only.



Western blot testing of human cerebral cortex lysate with APP antibody at 1ug/ml.  
Predicted molecular weight 79~120kDa depending on glycosylation level.

## Description

Additional name(s) for this target protein: Amyloid beta precursor protein; Abeta

## Application Notes

Optimal dilution of the APP antibody should be determined by the researcher.

## Immunogen

Amino acids DAEFRHDSGYEVHHQK were used as the immunogen for this APP antibody.

## **Storage**

Aliquot and store the APP antibody at -20oC.