



# **ARHGEF4 Antibody (R35175)**

Catalog No.	Formulation	Size
R35175-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	50649
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ARHGEF4 antibody is available for research use only.





Additional name(s) for this target protein: APC-stimulated guanine nucleotide exchange factor 4

# **Application Notes**

Optimal dilution of the ARHGEF4 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids PWEEPAGEKPSCS were used as the immunogen for this ARHGEF4 antibody.

## **Storage**

Aliquot and store the ARHGEF4 antibody at -20oC.