


## ARID5A Antibody (RQ6369)

Catalog No.	Formulation	Size
RQ6369	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

**Bulk quote request**

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit IgG
<b>Purity</b>	Antigen affinity purified
<b>Buffer</b>	Lyophilized from 1X PBS with 2% Trehalose
<b>UniProt</b>	Q03989
<b>Applications</b>	Western Blot : 0.5-1ug/ml Direct ELISA : 0.1-0.5ug/ml
<b>Limitations</b>	This ARID5A antibody is available for research use only.



Western blot testing of human 1) placenta and 2) K562 cell lysate with ARID5A antibody.  
Predicted molecular weight ~64 kDa.

## Description

Members of the ARID protein family, including ARID5A, have diverse functions but all appear to play important roles in development, tissue-specific gene expression, and regulation of cell growth.

## Application Notes

Optimal dilution of the ARID5A antibody should be determined by the researcher.

## Immunogen

Recombinant human protein (amino acids R55-L594) was used as the immunogen for the ARID5A antibody.

## Storage

After reconstitution, the ARID5A antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.