


B3GNT2 Antibody (R33827)

Catalog No.	Formulation	Size
R33827-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10678
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This B3GNT2 antibody is available for research use only.



Western blot testing of human pancreas lysate with B3GNT2 antibody at 0.1ug/ml. Predicted molecular weight ~46 kDa but may be observed at higher molecular weights due to glycosylation. Both observed bands are blocked by addition of immunizing peptide.

Description

Additional name(s) for this target protein: Beta-1,3-N-acetylglucosaminyltransferase 2

Application Notes

Optimal dilution of the B3GNT2 antibody should be determined by the researcher.

Immunogen

Amino acids EKHKGFRITFDIE were used as the immunogen for this B3GNT2 antibody.

Storage

Aliquot and store the B3GNT2 antibody at -20oC.