


# **BAIAP2 Antibody (R35908)**


Catalog No.	Formulation	Size
R35908-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10458
Localization	Cytoplasmic
Applications	Western Blot : 1-3ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This BAIAP2 antibody is available for research use only.



Western blot testing of human brain lysate with BAIAP2 antibody at 1 $\mu$ ml. Predicted molecular weight: ~57 kDa.



IHC testing of FFPE human brain (cortex) tissue with BAIAP2 antibody at 3.75ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

## **Description**

Additional name(s) for this target protein: BAI1-associated protein 2

## **Application Notes**

Optimal dilution of the BAIAP2 antibody should be determined by the researcher.

1. This BAIAP2 antibody is specific for isoform 1.

#### **Immunogen**

Amino acids SGSGTLVSTV were used as the immunogen for this BAIAP2 antibody.

#### **Storage**

Aliquot and store the BAIAP2 antibody at -20oC.