


BCL-2 Antibody (F47585)

Catalog No.	Formulation	Size
F47585-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47585-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine, Hamster, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P10415
Applications	Western Blot : 1:1000
Limitations	This BCL-2 antibody is available for research use only.



BCL-2 antibody pre-incubated without (Lane 1) and with (2) blocking peptide in 293 lysate

Description

This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Two transcript variants, produced by alternate splicing, differ in their C-terminal ends. [provided by RefSeq].

Application Notes

Titration of the BCL-2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 87-116 from the human protein was used as the immunogen for this BCL-2 antibody.

Storage

Aliquot the BCL-2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.