


BECN Antibody (F43508)

Catalog No.	Formulation	Size
F43508-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43508-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q14457
Applications	Western Blot : 1:1000
Limitations	This BECN antibody is available for research use only.



Description

Beclin 1 and its binding partner PI3K are required for the initiation of the formation of the autophagosome in autophagy.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the BECN antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 52-78 from the human protein was used as the immunogen for this BECN antibody.

Storage

Store at 4oC for up to one month. For long term, aliquot the BECN antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.