


# **BOC Antibody (R36257)**

Catalog No.	Formulation	Size
R36257-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	91653
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This BOC antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Brother of CDO

### **Application Notes**

Optimal dilution of the BOC antibody should be determined by the researcher.

#### **Immunogen**

Amino acids TQEVKTSGSSDRLR were used as the immunogen for this BOC antibody. **Storage** Aliquot and store the BOC antibody at -20oC.