


BRF2 Antibody (R35180)

Catalog No.	Formulation	Size
R35180-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55290
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This BRF2 antibody is available for research use only.



Western blot testing of HeLa lysate with BRF2 antibody at 1.5ug/ml. Predicted molecular weight: ~48 kDa.

Description

Additional name(s) for this target protein: BRF2 subunit of RNA polymerase III transcription initiation factor, BRF1-like

Application Notes

Optimal dilution of the BRF2 antibody should be determined by the researcher.

Immunogen

Amino acids QAARQAATSVPNPP were used as the immunogen for this BRF2 antibody.

Storage

Aliquot and store the BRF2 antibody at -20oC.