


BTBD7 Antibody (R33994)

Catalog No.	Formulation	Size
R33994-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55727
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This BTBD7 antibody is available for research use only.



Western blot of MDCK lysate overexpressing mouse Btbd7 (lane 3) tested with BTBD7 antibody. Lane 1: mock transfection, lane 2: full construct chemically suppressed

Description

Additional name(s) for this target protein: BTB/POZ domain-containing protein 7

Application Notes

Optimal dilution of the BTBD7 antibody should be determined by the researcher.

Immunogen

Amino acids ESKLYSLDHGHEKPQ were used as the immunogen for this BTBD7 antibody. **Storage** Aliquot and store the BTBD7 antibody at -20oC.