


# BTD Antibody (R33819)

Catalog No.	Formulation	Size
R33819-100UG	0.5~mg/ml in 1X TBS, pH7.3, with $0.5%$ BSA (US sourced) and $0.02%$ sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	686
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This BTD antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Biotinidase

### **Application Notes**

Optimal dilution of the BTD antibody should be determined by the researcher.

### **Immunogen**

Amino acids EKDAQEVHCDEATK were used as the immunogen for this BTD antibody.

# **Storage** Aliquot and store the BTD antibody at -20oC.