


C7orf28A Antibody (F40074)

Catalog No.	Formulation	Size
F40074-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40074-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Bovine
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P86791
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:50-1:100
Limitations	This C7orf28A antibody is available for research use only.

130 95 72 55 - ◀	C7orf28A antibody western blot analysis in mouse NIH3T3 lysate
36	
28	



C7orf28A antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human brain tissue.

Description

CG028 belongs to the UPF0550 family.

Application Notes

Titration of the C7orf28A antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 239-265 from the human protein was used as the immunogen for this C7orf28A antibody.

Storage

Aliquot the C7orf28A antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.