


CBF Antibody (F48151)

Catalog No.	Formulation	Size
F48151-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F48151-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q03701
Applications	Western Blot : 1:1000
Limitations	This CBF antibody is available for research use only.



Western blot analysis of CEBPZ / CBF antibody and NCI-H460 lysate. Predicted molecular weight: ~120kDa.

Description

CEBPZ stimulates transcription from the HSP70 promoter.

Application Notes

Titration of the CBF antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1001-1029 from the human protein was used as the immunogen for this CBF antibody.

Storage

Aliquot the CBF antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.