


CBF Antibody (F48151)

Catalog No.	Formulation	Size
F48151-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F48151-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q03701
Applications	Western Blot : 1:1000
Limitations	This CBF antibody is available for research use only.



Description

CEBPZ stimulates transcription from the HSP70 promoter.

Application Notes

Titration of the CBF antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1001-1029 from the human protein was used as the immunogen for this CBF antibody.

Storage

Aliquot the CBF antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.