


Cdk4 Antibody (F44483)

Catalog No.	Formulation	Size
F44483-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F44483-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P30285
Applications	Western Blot : 1:1000
Limitations	This Cdk4 antibody is available for research use only.



Description

Cdk4 is likely involved in the control of the cell cycle.

Application Notes

Titration of the Cdk4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 219-246 from the mouse protein was used as the immunogen for this Cdk4 antibody.

Storage

Aliquot the Cdk4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.