


Cebpd Antibody (F52321)

Catalog No.	Formulation	Size
F52321-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52321-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q03484
Applications	Western Blot : 1:1000
Limitations	This Cebpd antibody is available for research use only.



Cebpd antibody western blot analysis in 1) A549, 2) HeLa, 3) NCI-H460, 4) U-937, 5) mouse NIH3T3 cell line, rat 6) lung and 7) testis tissue lysate.

Description

CEBP is a DNA-binding protein that recognizes two different motifs: the CCAAT homology common to many promoters and the enhanced core homology common to many enhancers. Important transcriptional activator in the regulation of genes involved in immune and inflammatory responses, may play an important role in the regulation of the several genes associated with activation and/or differentiation of macrophages (By similarity). [UniProt]

Application Notes

Titration of the Cebpd antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 161-189 from rat CEBP delta was used as the immunogen for this Cebpd antibody.

Storage

Aliquot the Cebpd antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.