


# CHIT1 Antibody (R33796)

Catalog No.	Formulation	Size
R33796-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1118
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This CHIT1 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Chitinase 1

### **Application Notes**

Optimal dilution of the CHIT1 antibody should be determined by the researcher.

### **Immunogen**

Amino acids EHGPSPGQDTFCQGK were used as the immunogen for this CHIT1 antibody.

# **Storage** Aliquot and store the CHIT1 antibody at -20oC.