


Cntn4 Antibody (R35744)

Catalog No.	Formulation	Size
R35744-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	152330
Applications	Western Blot: 0.5-1.5ug/ml ELISA (peptide) LOD: 1:128000
Limitations	This Cntn4 antibody is available for research use only.



Western blot comparison of KO and WT tissue lysates with Cntn4 antibody at 0.5ug/ml. Lanes: KO mouse olfactory bulb (1) and cerebral cortext (2), and WT mouse olfactory bulb (3) and cerebral cortex (4). Lane 5: HEK293 lysate overexpressing mouse Cntn4.

Description

Additional name(s) for this target protein: Contactin 4

Application Notes

Optimal dilution of the Cntn4 antibody should be determined by the researcher.

Immunogen

Amino acids RRADGKPIARK were used as the immunogen for this Cntn4 antibody. **Storage** Aliquot and store the Cntn4 antibody at -20oC.