


# COL4A3BP Antibody (R33463)

Catalog No.	Formulation	Size
R33463-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

# **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10087
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This COL4A3BP antibody is available for research use only.



# **Description**

Additional name(s) for this target protein: Collagen, type IV, alpha 3 (Goodpasture antigen) binding protein

### **Application Notes**

Optimal dilution of the COL4A3BP antibody should be determined by the researcher.

# **Immunogen**

Amino acids DHDSAPLNNRC were used as the immunogen for this COL4A3BP antibody.

# **Storage** Aliquot and store the COL4A3BP antibody at -20oC.