


CT45A Antibody (F40015)

Catalog No.	Formulation	Size
F40015-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40015-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8NHU0
Applications	Western Blot : 1:1000
Limitations	This CT45A antibody is available for research use only.



Description

The function of this protein remains unknown.

Application Notes

Titration of the CT45A antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 100-129 from the human protein was used as the immunogen for this CT45A antibody.

Storage Aliquot the CT45A antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.