


## CT45A4 Antibody (F40019)

Catalog No.	Formulation	Size
F40019-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40019-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	P0DMU6
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This CT45A4 antibody is available for research use only.



CT45A antibody western blot analysis in WiDr lysate

## Description

Cancer/testis antigen family 45 member A4

## Application Notes

Titration of the CT45A4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 106-135 from the human protein was used as the immunogen for this CT45A4 antibody.

## **Storage**

Aliquot the CT45A4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.