


# **CT45A4 Antibody (F40019)**

Catalog No.	Formulation	Size
F40019-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40019-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P0DMU6
Applications	Western Blot: 1:1000
Limitations	This CT45A4 antibody is available for research use only.



### **Description**

Cancer/testis antigen family 45 member A4

## **Application Notes**

Titration of the CT45A4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### **Immunogen**

A portion of amino acids 106-135 from the human protein was used as the immunogen for this CT45A4 antibody.

# **Storage** Aliquot the CT45A4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.