


CTBP1 Antibody (F41046)

Catalog No.	Formulation	Size
F41046-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F41046-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q13363
Applications	Western Blot : 1:1000
Limitations	This CTBP1 antibody is available for research use only.



Western blot testing of CTBP1 antibody at 1:2000 dilution. Lane 1: HepG2 lysate; 2: SH-SY5Y lysate; 3: mouse brain lysate; 4: Raji lysate; Predicted molecular weight ~ 48 kDa.

Description

This gene encodes a protein that binds to the C-terminus of adenovirus E1A proteins. This phosphoprotein is a transcriptional repressor and may play a role during cellular proliferation. This protein and the product of a second closely related gene, CTBP2, can dimerize. Both proteins can also interact with a polycomb group protein complex which participates in regulation of gene expression during development. Alternative splicing of transcripts from this gene results in multiple transcript variants.

Application Notes

Titration of the CTBP1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 409-438 from the human protein was used as the immunogen for this CTBP1 antibody.

Storage

Aliquot the CTBP1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.