


# Cyt19 Antibody (R33106)

Catalog No.	Formulation	Size
R33106-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57412
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Cyt19 antibody is available for research use only.



IHC testing of FFPE human kidney with Cyt19 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## **Description**

Additional name(s) for this target protein: AS3MT; arsenic (+3 oxidation state) methyltransferase

### **Application Notes**

Optimal dilution of the Cyt19 antibody should be determined by the researcher.

### **Immunogen**

Amino acids AALRDAEIQKDVQ were used as the immunogen for this Cyt19 antibody.

# **Storage** Aliquot and store the Cyt19 antibody at -20oC.