


EH domain-containing protein 3 Antibody / EHD3 (F54929)


Catalog No.	Formulation	Size
F54929-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54929-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q9NZN3
Localization	Cell membrane, cytoplasmic
Applications	Flow Cytometry : 1:10-1:50 (1x10 ⁶ cells) Immunohistochemistry (FFPE) : 1:10-1:50 Western Blot : 1:500-1:1000
Limitations	This EH domain-containing protein 3 antibody is available for research use only.




IHC testing of FFPE human brain tissue with EH domain-containing protein 3 antibody.
HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.



Western blot testing of human brain tissue lysate with EH domain-containing protein 3 antibody. Predicted molecular weight ~61 kDa.



Western blot testing of mouse lung tissue lysate with EH domain-containing protein 3 antibody. Predicted molecular weight ~61 kDa.



Flow cytometry testing of human HepG2 cells with EH domain-containing protein 3 antibody; Blue=isotype control, Green= EH domain-containing protein 3 antibody.

Description

EHD3 plays a role in endocytic transport.

Application Notes

The stated application concentrations are suggested starting points. Titration of the EH domain-containing protein 3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 22-48 from the human protein was used as the immunogen for the EH domain-containing protein 3 antibody.

Storage

Aliquot the EH domain-containing protein 3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.

