

## EpCAM Antibody (R35320)

| Catalog No.  | Formulation   | Size   |
|--------------|---|--------|
| R35320-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |  |
|--------------------|--|
| Availability       | 1-3 business days  |
| Species Reactivity | Human  |
| Format             | Antigen affinity purified                                      |
| Clonality          | Polyclonal (goat origin)                                       |
| Isotype            | Goat Ig  |
| Purity             | Antigen affinity   |
| Gene ID            | 4072   |
| Applications       | Western Blot : 0.03-0.1ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations        | This EpCAM antibody is available for research use only.        |



Western blot testing of CACO-2 lysate with EpCAM antibody at 0.03ug/ml. Expected molecular weight: ~35 kDa (unmodified), 40-43 kDa (glycosylated).

## Description

Additional name(s) for this target protein: Epithelial cell adhesion molecule, TROP1

## Application Notes

Optimal dilution of the EpCAM antibody should be determined by the researcher.

## Immunogen

Amino acids QKEITTRYQLDPK were used as the immunogen for this EpCAM antibody.

## Storage

Aliquot and store the EpCAM antibody at -20oC.