


ETV5 Antibody (F47670)

Catalog No.	Formulation	Size
F47670-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47670-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P41161
Applications	Western Blot : 1:1000
Limitations	This ETV5 antibody is available for research use only.



Description

ETV5 binds to DNA sequences containing the consensus nucleotide core sequence GGAA.

Application Notes

Titration of the ETV5 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 8-37 from the human protein was used as the immunogen for this ETV5 antibody.

Storage

Aliquot the ETV5 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.