


## FAM159B Antibody (F53984)

| Catalog No.   | Formulation                                | Size    |
|---------------|--|---------|
| F53984-0.2ML  | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.2 ml  |
| F53984-0.05ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.05 ml |

[Bulk quote request](#)

|                           |   |
|---------------------------|---|
| <b>Availability</b>       | 1-3 business days   |
| <b>Species Reactivity</b> | Human   |
| <b>Format</b>             | Antigen affinity purified                                 |
| <b>Clonality</b>          | Polyclonal (rabbit origin)                                |
| <b>Isotype</b>            | Rabbit Ig   |
| <b>Purity</b>             | Antigen affinity  |
| <b>UniProt</b>            | A6NKG6  |
| <b>Applications</b>       | Western Blot : 1:1000-2000                                |
| <b>Limitations</b>        | This FAM159B antibody is available for research use only. |



Western blot testing of human 1) MKN45 and 2) PANC-1 cell lysate with FAM159B antibody at 1:2000. Predicted molecular weight: 18 kDa.

### Description

Additional name(s): Family with sequence similarity 159 member B

### Application Notes

Titration of the FAM159B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

A portion of amino acids 121-154 from the human protein was used as the immunogen for the FAM159B antibody.

## Storage

Aliquot the FAM159B antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.