


FAM159B Antibody (F53984)

Catalog No.	Formulation	Size
F53984-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53984-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	A6NKW6
Applications	Western Blot : 1:1000-2000
Limitations	This FAM159B antibody is available for research use only.



Western blot testing of human 1) MKN45 and 2) PANC-1 cell lysate with FAM159B antibody at 1:2000. Predicted molecular weight: 18 kDa.

Description

Additional name(s): Family with sequence similarity 159 member B

Application Notes

Titration of the FAM159B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 121-154 from the human protein was used as the immunogen for the FAM159B antibody.

Storage

Aliquot the FAM159B antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.