

FAM218A Antibody / C4orf39 (F54594)


Catalog No.	Formulation	Size
F54594-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54594-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity purified
UniProt	Q96MZ4
Applications	Western Blot : 1:500-1:2000 Immunofluorescence : 1:25 Flow Cytometry : 1:25 (1x10 ⁶ cells) Immunohistochemistry (FFPE) : 1:25
Limitations	This FAM218A antibody is available for research use only.




Western blot testing of human NCI-H460 cell lysate with FAM218A antibody. Predicted molecular weight ~17 kDa.



IHC testing of FFPE human pancreas tissue with FAM218A antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.



Immunofluorescent staining of human NCI-H460 cells with FAM218A antibody (green) and DAPI nuclear stain (blue).



Flow cytometry testing of human NCI-H460 cells with FAM218A antibody; Blue=isotype control, Green= FAM218A antibody.

Application Notes

The stated application concentrations are suggested starting points. Titration of the FAM218A antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 98-127 from the human protein was used as the immunogen for the FAM218A antibody.

Storage

Aliquot the FAM218A antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.