


FAM3A Antibody (F54169)

Catalog No.	Formulation	Size
F54169-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54169-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P98173
Applications	Western Blot : 1:1000-1:2000
Limitations	This FAM3A antibody is available for research use only.



Western blot testing of human testis tissue lysate with FAM3A antibody at 1:2000.
Predicted molecular weight ~25 kDa.

Description

Protein FAM3A may act as a defensin against invading fungal microorganisms.

Application Notes

The stated application concentrations are suggested starting points. Titration of the FAM3A antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 201-230 from human FAM3A was used as the immunogen for the FAM3A antibody.

Storage

Aliquot the FAM3A antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.