


FANCF Antibody (R35716)

Catalog No.	Formulation	Size
R35716-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2188
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This FANCF antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human FANCF-MYC probed with FANCF antibody [1ug/ml] in Lane A and anti-MYC [1/1000] in lane C. Mock-transfected HEK293 probed with FANCF antibody [1ug/ml] in Lane B. Predicted molecular weight: ~42kDa.

Description

Additional name(s) for this target protein: Fanconi anemia, complementation group F

Application Notes

Optimal dilution of the FANCF antibody should be determined by the researcher.

Immunogen

Amino acids RPQEELEPGIHK were used as the immunogen for this FANCF antibody.

Storage Aliquot and store the FANCF antibody at -20oC.						