


# FANCL Antibody (R35379)

Catalog No.	Formulation	Size
R35379-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55120
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This FANCL antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Fanconi anemia, complementation group L; FANCL

### **Application Notes**

Optimal dilution of the FANCL antibody should be determined by the researcher.

### **Immunogen**

Amino acids QNLKDVLEIDFP were used as the immunogen for this FANCL antibody. **Storage** Aliquot and store the FANCL antibody at -20oC.