


FBXO44 Antibody (R33107)

| Catalog No. | Formulation | Size |
|--------------|--|--------|
| R33107-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$ | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|--------------------|--|
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 93611 |
| Applications | Western Blot : 1.5-3ug/ml ELISA (peptide) LOD : 1:8000 |
| Limitations | This FBXO44 antibody is available for research use only. |



Description

Additional name(s) for this target protein: F-box protein 44

Application Notes

Optimal dilution of the FBXO44 antibody should be determined by the researcher.

Immunogen

Amino acids AALTPPEPPSAEP were used as the immunogen for this FBXO44 antibody. This sequence is specific to

Storage
Aliquot and store the FBXO44 antibody at -20oC.