


FGF10 Antibody (R30526)

Catalog No.	Formulation	Size
R30526	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	O15520
Applications	Western Blot : 0.5-1ug/ml
Limitations	This FGF10 antibody is available for research use only.



Western blot testing of FGF10 antibody and Lane 1: U87; 2: HeLa; 3: A549; 4: 293T; 5: HeLa cell lysate. Predicted molecular weight ~23 kDa.

Description

FGF 10, Fibroblast growth factor 10, is a protein that in humans is encoded by the FGF10 gene. The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF 10 contains 3 exons and spans at least 52.7 kb. By radioactive in situ hybridization, the FGF 10 gene is mapped to 5p13-p12. This protein exhibits mitogenic activity for keratinizing epidermal cells, but essentially no activity for fibroblasts. FGF 10 translocated into urothelial cell nuclei and initiated a signaling cascade that began with the heparin-dependent tyrosine phosphorylation of surface transmembrane receptors.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the FGF10 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

An amino acid sequence from the N-terminus of human FGF10 (DVRWRKLFSFTKYFLK) was used as the immunogen for this FGF10 antibody.

Storage

After reconstitution, the FGF10 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.