


FGF2 Antibody (F53229)

Catalog No.	Formulation	Size
F53229-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53229-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P09038
Applications	Western Blot : 1:2000
Limitations	This FGF2 antibody is available for research use only.



Western blot testing of FGF2 antibody at 1:2000 dilution; Lane 1: human A549 lysate, 2: mouse heart lysate; 3: rat lung lysate. Predicted molecular weight: 17-31 kDa (multiple isoforms).

Description

Plays an important role in the regulation of cell survival, cell division, angiogenesis, cell differentiation and cell migration. Functions as potent mitogen in vitro.

Application Notes

Titration of the FGF2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This FGF2 antibody was produced from a rabbit immunized with a recombinant protein.

Storage

Aliquot the FGF2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.